

Immobilizer (elbow)

☎ 00800 03032900

Brace to immobilize the elbow

Indications

- Post-surgical
- Paralysis
- Inflammatory rheumatic conditions

Features

- Simple alternative to wearing a cast
- Rigid shell made from PE with skinfriendly padding
- Broad velcro straps for fixation
- Physiological flexion angle of 58°

Optionally also available as elbow splint without padding and straps.

Elbow splint →
(without padding and straps)



Mode of Operation

The **elbow immobilizer** supports the treatment of the adult elbow and is used to immobilize the elbow.

Size chart and article numbers

Immobilizer

Size	Length of upper arm medial	Length of forearm medial	Art.-No.
S	18 cm	18 cm	37-2091
M	20 cm	20 cm	37-2092
L	22 cm	22 cm	37-2093

Elbow Splint

Size	Length of upper arm medial	Length of forearm medial	Art.-No.
S	18 cm	18 cm	37-2091S
M	20 cm	20 cm	37-2092S
L	22 cm	22 cm	37-2093S

Measurement chart on page 116



Measurement Chart

Noch nicht richtig übersetzt

