

Toe plate

☎ 00800 03032900


Product requiring individual adjustment by a technician

Indications

- Post-surgical fixation

Features

- Stable alternative to toe slings
- Rigid shell made out of PE
- Snuggly velours strap
- Individual fixation



Individual cord set
toe plate

Mode of Operation

The **Toe plate** assists treatment in the application area of the foot. It is used to immobilize the toes.

Size chart and article numbers

| Nmae | Art.-No. |
|-------------------------------|-----------|
| Toe plate with cord | 37-8000 |
| Toe plate with strap | 37-8500 |
| Individual cord set toe plate | 770-026-S |

Measurement Chart

Noch nicht richtig übersetzt

