

Immobilizer (foot)

Brace to immobilize the foot

Indications

- Post-surgical to immobilize
- Post-surgical after achill tendon surgery
- Shortening of the plantar fascia
- Prophylaxis against talipedes equinus
- Pain from overstrain

Features

- Simple alternative to wearing a cast
- Rigid shell made from PE with skinfriendly padding
- Broad velcro straps for fixation
- Adjustable angle

Optionally also available as foot splint without padding and straps.

Foot splint →
(without padding
and straps)



Mode of Operation

The **foot immobilizer** supports the treatment of the adult foot and is used to immobilize the foot.

Size chart and article numbers

Immobilizer

Size	EU-shoe size	Sole Length	Art.-No.
S	36 - 39	25 cm	37-4051
M	40 - 43	27 cm	37-4052
L	44 - 47	30 cm	37-4053

Foot Splint

Size	EU-shoe size	Sole Length	Art.-No.
S	36 - 39	25 cm	37-4051S
M	40 - 43	27 cm	37-4052S
L	44 - 47	30 cm	37-4053S



Measurement Chart

Noch nicht richtig übersetzt

