

Finger Brace

Patents EP 3 421 016

Indications

- Joint contractures due to
 - burns
 - stroke
 - spinal injury
 - paralysis
 - craniocerebral injury
 - surgery
 - arthrosis and chronic polyarthritis
- cut injury
- tendon suture
- Morbus Dupuytren
- to prevent recurring contractures following arthrolysis

Features

- Extension and flexion therapy
- Individually adjustable spring force
- Individually adjustable shell and strap system
- Easy handling and comfortable in use
- Adjustable hand plate
- No tool necessary for adjusting spring force
- Individual toolless adaptation of chains to finger length



Dynamic
splints

Mode of Operation

The **Manusan** is used to treat an extension deficit of the finger joints. Active use trains the appendicular finger muscles. Passive use stimulates growth of contracted tissue due to continuous dynamic tension.

Size chart and article numbers

Name	Circumference	Width Palm plate	Art.-No. left	Art.-No. right
Manusan L/L	22 - 30 cm	10 cm	955LL-L	955LL-R
Manusan M/M	18 - 24 cm	9 cm	955MM-L	955MM-R
Thumb element				956

Other sizes available on request.

Measurement Chart

Noch nicht richtig übersetzt

