

# MKS® Core Compression Bandage

## Compression bandage

### Indications

- Prevention of blood congestion
- Support during transfer

### Features

- Semi-elastic bandaging material
- Abdominal piece with hook and loop fastener
- Individual strap course
- Detachable lateral double straps

### Mode of Operation

The function of the MKS® Core Compression Bandage is based on the principle of compression. The stomach part, which consists of several straps can be individually adapted to the patient's needs by means of a side binding piece. This produces a compression pressure on the lower abdomen and can thus counteract blood retention.



### Size chart and article numbers

#### MKS® Core Compression Bandage (4 straps)

Size	Body circumference	Abdominal height	Height of back piece	Art.-No.
1	60 - 70 cm	32 cm	32 cm	583-001-S
2	70 - 80 cm	32 cm	32 cm	583-002-S
3	80 - 90 cm	34 cm	34 cm	583-003-S
4	90 - 105 cm	34 cm	34 cm	583-004-S
5	105 - 120 cm	34 cm	34 cm	583-005-S
6	120 - 135 cm	34 cm	34 cm	583-006-S

#### MKS® Core Compression Bandage (3 straps)

Size	Body circumference	Abdominal height	Height of back piece	Art.-No.
1	60 - 70 cm	22 cm	22 cm	583-001N-S
2	70 - 80 cm	22 cm	22 cm	583-002N-S
3	80 - 90 cm	22 cm	22 cm	583-003N-S
4	90 - 105 cm	22 cm	22 cm	583-004N-S
5	105 - 120 cm	22 cm	22 cm	583-005N-S
6	120 - 135 cm	22 cm	22 cm	583-006N-S

Optionally side binding pieces for additional compression.



Optionally side binding pieces

