

Lumbar bandage

Indications

- Lower back pain
- Moderately severe lumbar sciatica with slight failures in intervertebral disc protrusion/prolapse
- Moderately severe radicular lumbar syndrome, pseudo-radicular lumbar syndrome
- Spondylolisthesis, grade I with lumbago
- Moderately severe lumbar deformity (relapsing treatment, complicated course)
- Condition after intervertebral disc surgery with small up to medium-sized prolapse
- Spondylarthrosis
- Sacroiliac joint irritation
- Lumbar muscle insufficiency

Features

- Semi-elastic bandaging material
- Lateral reinforcing stays
- Firm abdominal piece with hook and loop fastener
- Bridging rods (spring band steel)
- Stable bridging frame that adjusts to contour, with window



Mode of Operation

The MKS® Lumbo fenestra functions according to the bridging principle. The wide bandage with the additional lateral double straps straightens and hence relieves the lumbar vertebral column. The stable bridging frame that adjusts to contour contains a window for easy access to the wound and pressure relief of the injured area.

Size chart and article numbers

Size	Body circumference	Abdominal height	Height of back piece	Art.-No.
1	60 - 70 cm	18 cm	32 cm	790-11
2	70 - 80 cm	18 cm	32 cm	790-21
3	80 - 90 cm	18 cm	34 cm	790-31
4	90 - 105 cm	18 cm	34 cm	790-41
5	105 - 120 cm	18 cm	38 cm	790-51
6	120 - 135 cm	18 cm	38 cm	790-61



