

# taluai<sup>®</sup> Ankle Brace

## Ankle orthosis for stabilisation in one plane

### Indications

- Post-surgical protection after ligament repair/reconstruction
- Conservative therapy of ankle joint distortions of all degrees of severity
- Functional therapy of ligament ruptures
- Ankle joint fractures (Weber A)
- Ligament instabilities
- Ligament injuries of the ankle joint
- Severe ankle joint distortion with loosened ligament
- Conservative or functional follow-up treatment of ankle joint fractures and ligament ruptures

### Features

- Air chambers made out of skin-friendly polyurethane sheeting
- PVC-free and free of plasticizers
- Air chambers inner core made out of PU foam filling
- Available in white and black

**Optional** pronation strap for upgrading to talupro<sup>®</sup>

### Mode of Operation

The taluai<sup>®</sup> provides medial and lateral support to the ankle joint while allowing a great degree of flexion and extension. This allows for a normal motion sequence. The air chambers cushion the supports to prevent unpleasant pressure. Stronger padding at both ankles provides higher pressure at the bottom third of the supporting area. The applied pressure decreases from the bottom towards the top, helping to reduce swelling more rapidly.

### Size chart and article numbers

Size	Body height	Colour	Art.-No. left	Art.-No. right
<b>S</b>	under 158 cm	white	600S-L	600S-R
<b>M</b>	158 - 163 cm	white	600M-L	600M-R
<b>L</b>	from 163 cm	white	600L-L	600L-R
<b>M</b>	to 163 cm	black	600MS-L	600MS-R
<b>L</b>	from 163 cm	black	600LS-L	600LS-R
<b>Pronation strap for upgrading to talupro<sup>®</sup></b>		white	481-020-S	481-025-S
<b>Pronation strap for upgrading to talupro<sup>®</sup></b>		black	482-020-S	482-025-S



# Measurement Chart

Noch nicht richtig übersetzt

